

APPLICATION MANUAL

1. Kit Content

Cat no	Content
ASFV E-55	1. PCR-reagent 1 tube with black cap (includes enzyme in buffer) – 0.55 ml
	2. ASFV Primers 1 tube with yellow cap (includes primers and probe) – 0.55 ml
	3. Reagent “DNA-express” 1 tube with violet cap (includes special sorbent in buffer) – 1.37 ml
	4. External control (EC) 1 tube with green cap – 0.55 ml
	5. Positive control (PC) 1 tube with red cap – 0.5 ml
	6. Negative control (NC) 2 tubes with white caps – 1 ml

2. Intended Use

ASFV Kit is used for the detection of African Swine Fever Virus (ASFV) DNA in serum, plasma and tissue samples from pig and wild boar by using real-time PCR instruments.

3. Storage

The components of the -ASFV PCR Kit could be stored at +2 - +4°C for up to 12 months from the date of manufacturing without reduction in sensitivity.

4. Introduction

African swine fever is an infectious disease of domestic and wild pigs of all breeds and ages, caused by a virus (member of the Asfarviridae family, genus Asfivirus) that produces a range of syndromes, varying from peracute, acute to chronic disease and apparently healthy virus carriers. Pigs are the only domestic animal species that is naturally infected by ASFV. European wild boars and feral pigs are also susceptible to the disease, exhibiting clinical signs and mortality rates similar to those observed in domestic pigs.

The -ASFV PCR Kit is a highly sensitive and specific solution for the detection of ASFV DNA in samples from pigs and wild boars. The express method of DNA purification from the wide range of clinical samples allows detecting the pathogen early and easily.

5. Principle

Polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time PCR, the amplified product is detected using fluorescent dyes. Monitoring the fluorescence intensities during the PCR run (i.e., in real time) allows detection of the accumulating product without the need to re-open the reaction tubes afterward. The -ASFV PCR Kit contains all of the necessary reagents for the detection of ASFV DNA, including reagents for Express DNA purification, external, positive and negative controls. An external control excludes the possibility of false-negative results.

The kit uses three specific primer/probe combinations: Two - for ASFV DNA yielding HEXTM fluorescence and one for an external control yielding FAMTM fluorescence. A PC serves to verify the functionality of the reaction mixture for the amplification of the ASFV DNA target.

6. Equipment and reagents to be supplied by user

- Pipets (0.5 µl – 1000 µl);
- Nuclease-free aerosol-resistant pipet tips with filters;
- Sterile 1.5 ml tubes;
- Benchtop centrifuge with rotor for 1.5 ml tubes (RCF \geq 10,000 × g);
- Vortex mixer;
- Disposable gloves, powderless;
- Biohazard waste container;
- Refrigerator and freezer;
- Tube racks;
- Real-time PCR machine (e.g. Rotor-Gene 3000, 6000, Q; Biorad CFX 96; Agilent MX 3005P and etc.);
- PCR tubes, suitable for used PCR machine.

While working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles.

7. Sample preparation

-Plasma and serum samples are used without special preparation.

-Tissue and organ samples should be prepared as 10% suspension (5-20% range is allowed) in phosphate buffer saline (PBS) and then homogenized using porcelain mortars and pestles or automatic homogenizers (recommended). After homogenization the suspension should be transferred to a 1.5 ml tube and centrifuged at $10,000 \times g$ for 5 min. The supernatant is used for DNA extraction with reagent "DNA-Express".

8. Protocol for DNA purification using the reagent "DNA-express"

1. Prepare 1.5 ml or 2 ml microcentrifuge tubes in count of tested samples but adding one extra tubes for NC extraction. Mark the tubes carefully.

2. Thoroughly mix the reagent "DNA-express" by vortexing before use.

3. Pipet 25 μ l of the reagent "DNA-express" into each tube. Then add 5 μ l of EC into each one.

3. Add 25 μ l of the test material (plasma, serum, tissue extract) into the corresponding tube.

4. Add 25 μ l of NC into a separate tube and mark it as Negative Control of Extraction (NCE).

5. Close the reaction tubes with the corresponding caps. Pulse-vortex tubes for approximately 5 seconds.

6. Incubate at 95°C for 10 min with constant mixing (more than 1000 rpm recommended).

7. Centrifuge for 5 min at minimum $10,000 \times g$.

8. Carefully transfer the supernatant (about 30 μ l) into the marked 1.5 ml microcentrifuge tubes. Store the extracted DNA at -20°C until required for PCR amplification procedure.

9. Protocol for real-time PCR detection of ASFV

1. Before use mix reagents by inverting tubes several times and spin the reagents briefly.

2. In a 1.5 ml microcentrifuge tube prepare the PCR reaction "Master Mix" in bulk for the number of samples to be assayed but allowing for one extra sample:

(10 μ l PCR-reagent + 10 μ l ASFV Primers)*(N+1),

where N is number of samples.

3. Pipet 20 μ l of the "Master Mix" into each reaction tube. Then add 10 μ l of the sample DNA (Table 1). Include positive and negative control reactions.

Positive Control: Use 10 μ l of the PC instead of sample DNA.

Negative Control: Use 10 μ l of the NC instead of sample DNA.

Table 1. Preparation of reaction mix

Component	Volume
Master Mix	20
Sample	10
Total volume	30

- Close the reaction tubes with the corresponding caps.
- Set the filters for the reporter dyes in the software of your thermal cycler according to Table 2.

Table 2. Filter settings for the reporter

ASFV/External control	Reporter
ASFV	HEX/Yellow
External Control	FAM/Green

- Run the real-time PCR protocol according to Table 3.

Table 3. PCR protocol

Temperature	Time	Number of cycles
95°C	15 min	1
95°C	5 s	10
60°C	10 s	35
67	10 s	
95°C	5 s	
60°C	10 s, Fluorescence data collection (FAM/GREEN; HEX/YELLOW)	
67	10 s	

10. Interpretation of Results

For the assay to be valid, the PC must give a signal in both the FAM and HEX channels with a $Ct \leq 27$. The NC and NCE must give no signal in HEX channel.

The possible results of control samples are also summarized in Table 4.

Table 4. Evaluation of control samples

Controls	FAM	HEX
NCE	$Ct \leq 27$	No Signal
NC	No Signal	No Signal
PC	$Ct \leq 27$	$Ct \leq 27$

The following results are possible if working with unknown samples (Table 5).

The sample is positive for ASFV, and the assay is valid, if the following criteria are met:

- The sample yield a signal in HEX channel;
- The NC and NCE does not yield a signal in the HEX channel.

The sample is negative for ASFV, and the assay is valid, if the following criteria are met:

- The sample yields a signal in the FAM channel but not in the HEX channel;
- The NC and NCE does not yield a signal in the HEX channel.

The sample results are inconclusive, and the assay is invalid, if the following occurs:

- The sample yields no signal in the FAM and HEX channels.

Table 5. Results interpretation table.

Sample result	FAM (EC)	HEX (pathogen)
ASFV positive	+/-	+
ASFV negative	+	-
Inconclusive	-	-

11. Troubleshooting

Trouble	Possible causes	Solution
No signal for the EC	PCR inhibition	-dilute the DNA 1:5 in nuclease free water; -repeat the whole test procedure starting with new sample material; -use another DNA extraction method
No signal for the PC	incorrect preparation of the reaction mixture	carefully prepare new reaction mixture
	incorrect cycling conditions	carefully set cycling conditions
	improper storage	use properly stored kit
NC signal presence	contamination	-perform decontamination procedures; - use new kit

AUTHORIZED BY
YZZY Biotech
R&D laboratory