

### Components:

-YZZY Start DNA Polymerase 5U/ul (5000 U)  
10x1 ml

### Introduction

YZZY Start is a thermostable, recombinant polymerase with aptamer mediated, reversible hot start & immediate activation.

Similar to antibody-based methods, the enzyme is activated by heating. Dropping the temperature below +39°C acts like a switch, shutting off polymerase activity. At temperatures above +45°C, YZZY Start Polymerase is fully activated. ideal for specific priming and fast PCR.

The enzyme catalyzes the template-dependent polymerization of nucleotides into duplex DNA in the 5'→3' direction and also exhibits 5'- exonuclease activity.

- High-specificity PCR;
- Routine PCR;
- Real-Time PCR with dual labeled probes and intercalating dyes.

### Protocol

#### Important point before starting:

To reduce the pipetting error when setting several parallel reactions, it is recommended to prepare a common PCR premix containing all the components of the mixture (water, dNTP, buffer, DNA polymerase, primers) with the exception of the DNA template. The volume of PCR premix is calculated based on the number of reactions with the addition of one additional sample.

1. Thaw all reagents, template DNA and primers
2. Prepare a reaction mix according to Table1 (For routine PCR) and Table2 (For Real-Time PCR with Taqman probe).

**Table 1. Reaction setup for routine PCR**

Component	Volume/ 25 ul reaction	Final concentration
<b>10X Buffer<sup>1</sup></b>	2,5 ul	1X
<b>Forward primer (10 uM)</b>	0,5 - 1 ul	0,2 - 0,4 uM
<b>Reverse primer (10 mM)</b>	0,5 - 1 ul	0,2 - 0,4 uM
<b>dNTPs (10 mM)</b>	0,5 ul	200 uM
<b>MgCl<sub>2</sub> (50 mM)<sup>2</sup></b>	1 – 1,25 ul	2 – 2,5 mM
<b>Template DNA</b>	Variable	0,1 ng - 250 ng (for high complexity genomic DNA)

<b>YZZY Start U/ul</b>	0,1 – 0,2 ul	0,5 - 1 U
<b>Nuclease Free H2O</b>	to 25 ul	

1 - Any standard PCR buffer can be used.

As a starting buffer for optimization, we recommend using a buffer of the following composition 10X AR: 650 mM Tris-HCl (pH 8.8), 166 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 0.5% Tween 20.

**Table 2. Reaction setup for PCR with Taqman probe**

Component	Volume/ 25 ul reaction	Final concentration
<b>10X Buffer<sup>1</sup></b>	2,5 ul	1X
<b>Forward primer (10 uM)</b>	0,7 - 1 ul	0,28 - 0,4 uM
<b>Reverse primer (10 uM)</b>	0,7 - 1 ul	0,28 - 0,4 uM
<b>Probe (10 uM)</b>	0,25 – 0,7 ul	0,1 – 0,28 uM
<b>dNTPs (10mM)</b>	0,5 ul	200 uM
<b>MgCl<sub>2</sub> (50mM)<sup>2</sup></b>	2 ul	4 mM
<b>Template DNA</b>	Variable	0,1 ng - 250 ng (for high complexity genomic DNA)
<b>YZZY Start 5 U/ul</b>	0,4 ul	2 U
<b>Nuclease Free H2O</b>	to 25 ul	

1 - As a starting buffer for optimization, we recommend using a buffer of the following composition **10X A**: 650 mM Tris-HCl (pH 8.8), 200 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 0.5% Tween 20.

- Mix the reaction thoroughly, and dispense appropriate volumes into PCR tubes or plates.
- Add template DNA to the individual PCR tubes or wells containing the reaction mix.
- Program your thermocycler according to the program outlined in Table 3 (for routine PCR) and Table 4 (For Real-Time PCR).

**Table 3. PCR conditions for routine PCR**

Cycle Step	Tempera-ture	Time	Num-ber of cycles
<b>Initial Denaturation</b>	95 °C	2 min	1
<b>Denaturation</b>	95 °C	15-30 s	25-35
<b>Annealing</b>	55°C – 67 °C	15-30 s	
<b>Elongation</b>	67 °C	1min/kb	
<b>Final extension</b>	67 °C	5min	1

The recommended extension temperature is 67°C. Extension times are generally 1 minute per kb. A final extension of 5 minutes at 67°C is recommended.

Generally, 25–35 cycles yield sufficient product. Up to 45 cycles may be required to detect low-copy-number targets.



**Table 4 Real-time cycler conditions**

Cycle Step	Tempera-ture	Time	Num-ber of cycles
<b>Initial Denaturation</b>	95 °C	2 min	1
<b>Denaturation</b>	95 °C	10 s	45
<b>Annealing</b>	60 °C	15 s + Plate Read	
<b>Elongation</b>	67 °C	15 s	

### Shipping and storage

**YZZY Start DNA polymerase** is shipped at ambient temperature (below 25°C) up to 30 days, at +4°C for a long time.

**YZZY Start DNA polymerase** upon arrival should be stored at –20°C, and has a shelflife of 24 months when stored properly under these conditions.

**YZZY Start polymerase** in the supplied buffer can be stored at room temperature (25 ° C) for 1 month without significant changes in activity and inactivation level.

PREPARED BY „YZZY BIOTECH“  
R&D LABORATORY

